

What's Your Handle?



0010



AB03



AB04



0012

0011

0001

AB01

AB05

Aaron Medical manufactures a diverse line of ophthalmic specialty products. Among these are three ophthalmic burr power handles and five different replacement burr tips.



ISO 9001 CERTIFIED
AARON
MEDICAL INDUSTRIES
a *Bovie* company

7100 30th Ave. N. • St. Petersburg, FL 33710-2902
Phone 800-537-2790 • 727-384-2323
Fax 727-347-9144 • www.aaronmed.com

20M 3/00 MC17354 Rev.0


Aaron Burr Handles & Tips Descriptions



Aaron Burr 2,

Catalog reference number 0010

This Aaron burr 2 is a deluxe power handle, and is designed for corneal rust ring removal, polishing the pterygium bed after surgical removal, and for lid margin lesions. This unit comes with one 0001 replacement tip, AA alkaline batteries and a hard plastic case. There are three placement tips available for this unit the 0001, 0011, and 0012.



0001



0011



0012

Replacement Tip

Catalog reference number 0001

The 0001 is a 1mm replacement tip and when used with the Aaron Burr two is ideal for corneal rust ring removal. These tips are disposable and packaged individually sterile and sold 10 per box.

Replacement Tips

Catalog reference numbers 0011 & 0012

The 0011 is a 2.5mm diamond burr tip and the 0012 is a 5mm diamond burr tip. When used with the Aaron Burr 2 (#0010) they may be used for polishing the pterygium bed after surgical removal and for lid margin lesions. These tips are autoclavable and are sold individually.

Aaron Burr 3,

Catalog reference number AB03

This power handle is our most economical model and is designed strictly for corneal rust ring removal. This handle is made of durable plastic and has a push button activation switch.



Aaron Burr 4,

Catalog reference number AB04

The Aaron Burr 4 is also designed for corneal rust ring procedures is chrome and also comes with a push button activation switch.



AB01



AB05

Replacement Tips

Catalog reference number AB01 & AB05

Both the AB03 and AB04 have two tip replacement possibilities. These are the AB01, 1mm replacement tip and the AB05, 1/2mm replacement tip. These tips are disposable and packaged individually sterile and sold 10 per box.

Distributed By: