


# DISPOSABLE ELECTRODES FROM A NAME YOU TRUST!



All Aaron® blade, ball, needle, loop, and square disposable electrodes use the standard 3/32" diameter stainless steel shafts. These electrodes, like all Aaron® electrodes, feature safety grip insulators, combining patient and user safety with easy insertion into and removal from the surgical pencil.

The loop and square electrodes feature tungsten wire for superior shape and wire integrity throughout the excision procedure, while minimizing thermal damage and yielding excellent histopathology specimens. Each molded loop and square electrode includes our patented tray to protect the loops from bending during shipping and handling.

Ball electrodes in 3mm and 5mm sizes are available for fulguration and desiccation, while needle electrodes are available for precise, pinpoint cutting and coagulation.

SuperCut™ Needles (ES 60-63) are made of Tungsten and keep their shape throughout the procedure.

The Aaron Arthroscopic electrodes are available in four tip styles and designed to fit all standard electro-surgical pencils.



7100 30th Avenue North  
 St. Petersburg, FL 33710-2902  
 U.S. Ph 1-800-537-2790 • Fx 1-800-323-1640  
 Int'l. Ph +1-727-384-2323 • Fx +1-727-344-3876  
[www.aaronmed.com](http://www.aaronmed.com) • [sales@aaronmed.com](mailto:sales@aaronmed.com)

# LIST OF STYLES AND CONVERSION CHART

## Loop Size Locator

	STRYKER/COLORADO	PREMIER	WALLACH	CONMED	VALLEY	ACMI/CI/RCON	UTAH		
SuperCur™ Tungsten Needles	ES60	3CM Super Fine Needle, 5/box	N101A						
	ES61	3CM Modified Super Fine Needle, 5/box	N103A		E1651				
	ES62	4.5CM Super Fine Needle, 5/box	N102A						
Needle Electrodes	ES63	4.5CM Modified Super Fine Needle, 5/box	N104A						
	ES02	Standard Needle, 25/box	COOPER	PREMIER	WALLACH	CONMED	VALLEY	ACMI/CI/RCON	UTAH
	ES03	Extended Needle, 25/box			909135	275-1502/9113	E1552		
Dermal Loops	ES38	Modified Needle, 25/box				275-1504/9116	E1552-6		
	ES40	Extended Modified Needle, 25/box	N2512			138025		6162901	
	ES22	1/4" Short Shaft Loop, 5/box							
	ES23	5/16" Short Shaft Loop, 5/box							
Dermal Tips	ES24	3/8" Short Shaft Loop, 5/box							
	ES25	1/2" Short Shaft Loop, 5/box							
	ES26	3/4" Short Shaft Loop, 5/box							
	A804	Sharp Dermal Tip, Non-Sterile, 100/box							
Blade Electrodes	A805	Sharp Dermal Tip, Sterile, 50/box							
	A806	Blunt Dermal Tip, Non-Sterile, 100/box							
	A807	Blunt Dermal Tip, Sterile, 50/box							
	ES01	Standard Blade, 50/box					275-1500/9111	E1551G	
Ball Electrodes	ES04	Extended Blade, 25/box					275-1503/9114	E1551-6	
	ES37	Modified Blade, 25/box					275-1524	E1553	
	ES39	Extended Modified Blade, 25/box					275-1525		
	ES18	Anoled Blade, 5/box							
LLETZ Loop Electrodes	ES06	Extended 3mm Ball, 5/box	B0312			60-5181-102	E1563	6156901	DBL-311
	ES07	Extended 5mm Ball, 5/box	B0512			60-5181-103	E1564	6157901	DBL-511
	ES20	3/16" Ball, 25/box				275-1501/9112	E1550		
	ES21	5/32" Ball, 25/box							
Arthroscopic Electrodes	ES50	1/8" Ball, 25/box							
	ES08	5mm x 5mm Loop, 5/box						6133901	DLP-T11
	ES09	10mm x 10mm Loop, 5/box	R1010					6134901	DLP-S11
	ES10	15mm x 8mm Loop, 5/box	R1507*					6135901	
	ES11	15mm x 10mm Loop, 5/box				60-5180-102		6136901	DLP-M11
	ES12	20mm x 8mm Loop, 5/box	R2008					6137901	DLP-L11
	ES13	20mm x 15mm Loop, 5/box	R2015					6138901	
	ES31	20mm x 20mm Loop, 5/box	R2020						
	ES42	20mm x 12mm Loop, 5/box							
	ES43	15mm x 15mm Loop, 5/box							
	ES44	15mm x 5mm Loop, 5/box	R1515						
	ES44	15mm x 5mm Loop, 5/box	R1505						
	ES45	13mm x 13mm Loop, 5/box							
	ES46	10mm x 5mm Loop, 5/box							
	ES47	25mm x 10mm Loop, 5/box							
	ES48	10mm x 3mm Loop, 5/box	R2510						
ES49	20mm x 10mm Loop, 5/box								
ES51	13mm x 8mm Loop, 5/box	R2010							
ES52	10mm x 8mm Loop, 5/box								
ES52	10mm x 8mm Loop, 5/box					60-5180-104			
ES52	15mm x 12mm Loop, 5/box					60-5180-101*			
ES53	15mm x 8mm Loop, 5/box								
ES57	25mm x 8mm Loop, 5/box								
ES57	25mm x 8mm Loop, 5/box								
ES58	25mm x 22mm Loop, 5/box								
ES14	5mm x 5mm Square, 5/box							6153901	
ES15	10mm x 4mm Square, 5/box	S1004						6154901	DLP-SQ2
ES16	10mm x 8mm Square, 5/box							6155901	DLP-SQ1
ES41	5mm x 10mm Square, 5/box								
ES56	10mm x 10mm Square, 5/box								
AR00	Electrode Hook 90, 5/box								E1510
AR01	Electrode Hook 45, 5/box								E1511
AR02	Electrode Acromioplasty, 5/box								E1512
AR03	Electrode Meniscectomy, 5/box								E1513

\* shows 1mm difference in size

