

Ophthalmic Products Made With Vision!



Aaron manufactures seven different cauteries specifically designed for ophthalmology. These patented Aaron RAM™ cauteries have an impressive four-year shelf life and are quite simply the best cauteries ever made! With a variety of low temperatures and fine tips, the surgeon is assured of precise pinpoint hemostasis and minimal tissue damage.

Additionally, this variety allows the surgeon to pick the right cautery at the right temperature, with confidence. All disposable cauteries come individually packaged sterile, 10 per box.



Aaron Medical manufactures a diverse line of ophthalmic specialty products. Five different specialty lights, for a variety of applications, are available. The cobalt blue spot light (6802A) and slit image lamp (6803A) are for use with fluorescein dye. Also available is a white slit image light (6801A). The two other specialty lights are the Tech-Lite (6800A), which is an oversized penlight and the (6804A) used specifically for eye muscle tests. All of the above specialty lights have a four-hour burn time and a one-hand activating switch.

The (0002) sterile eye bubble consists of a clear protective dome window with an adhesive perimeter for easy attachment. The bubble keeps the cornea hydrated and is also excellent protection for the eye during surgery.



7100 30th Avenue North
St. Petersburg, FL 33710-2902
800-537-2790 • Fax 800-323-1640
(727) 384-2323 • Fax (727) 347-9144
www.aaronmed.com
E-mail=customer.service@aaronmed.com

Aaron Medical's Line of Ophthalmic Products

CATALOG NUMBER

DESCRIPTION


- AA00 Fine Tip Low Temp
- AA02 Fine Tip Low Temp Adjustable
- AA04 Elongated Fine Tip Low Temp
- AA90 Fine Tip Micro Low Temp

- HIT0 Replaceable Battery and Tip Low Temp
- H100 Low Temp Fine Replacement Tip with Sterile Drape
- H104 Low Temp Elongated Fine Replacement Tip with Sterile Drape

Ophthalmic Cauteries

Aaron manufactures six different disposable cauteries designed especially for ophthalmology. These low temperature cauteries have the proven "USA made" design with a low temperature for precise pinpoint hemostasis and minimum tissue damage. By having five ophthalmic cauteries, Aaron allows the surgeon to pick with confidence the right cautery at the right temperature. Packaged individually sterile, 10 per shelf pack.

Also available is our low temperature replaceable battery and tip cautery.



- 0002 Sterile Eye Bubble

Aaron Eye Bubble

This sterile, self-adhesive dressing consists of a clear dome window, allowing the patient unobstructed vision. The Eye Bubble is available in a comfortable, protruding style. In most cases, the physician can perform an eye examination without removing the bandage. In conditions where drying is a factor, the bandage effectively keeps the cornea hydrated and helps protect the eye during anesthesia and during surgery. Individually packaged sterile, 10 per box.



- 6800A Tech Light
- 6801A White Slit
- 6802A Cobalt Blue
- 6803A Cobalt Slit
- 6804A 1mm Muscle Test

Ophthalmic Lighting

Each ophthalmic light contains a constant on/off switch, two AA alkaline batteries, a four hour burn time and a four year shelf life. The lights come in the following styles: Cobalt Blue, Cobalt Slit, White Slit, Standard, and 1mm for eye muscle tests. Individually packaged.

- 6303 Cobalt Blue Penlight 3 Pack
- 6666 Disposable Penlight 6 Pack
- 666R Replaceable Battery and Lamp Penlight 6 Pack

Penlights

Disposable Penlights have a 1 year shelf life and are available with a cobalt blue filter for ophthalmic procedures. Also available in replaceable battery and lamp format.

Distributed by: