











**MEG™ Handle with ratchet  
Insulated with power post  
Reusable MIHNDI**

-  Aggressive Grasper  
Disposable Cartridge MI32AGG-DC  
Reusable Tip Insert MI32AGG-RI
-  Babcock Grasper  
Disposable Cartridge MI32BCG-DC  
Reusable Tip Insert MI32BCG-RI
-  Atraumatic Grasper  
Disposable Cartridge MI32ATG-DC  
Reusable Tip Insert MI32ATG-RI
-  Allis Grasper  
Disposable Cartridge MI32ALG-DC  
Reusable Tip Insert MI32ALG-RI
-  Curved Metzenbaum Scissors  
Disposable Cartridge MI32CMS-DC  
Disposable Tip Insert MI32CMS-DI

-  Blunt Dissector  
Disposable Cartridge MI32BLD-DC  
Reusable Tip Insert MI32BLD-RI
-  Dolphin Nose Dissector  
Disposable Cartridge MI32DND-DC  
Reusable Tip Insert MI32DND-RI
-  Maryland Dissector  
Disposable Cartridge MI32MYD-DC  
Reusable Tip Insert MI32MYD-RI
-  Mixer Dissector  
Disposable Cartridge MI32MXD-DC  
Reusable Tip Insert MI32MXD-RI
-  Straight Scissors  
Disposable Cartridge MI32STS-DC  
Disposable Tip Insert MI32STS-DI

*For more information concerning MEG™ line  
of Advanced Endoscopic Instruments,  
please contact your Local Distributor.*



**BOVIE MEDICAL CORPORATION**  
7100 30<sup>th</sup> Avenue North • St. Petersburg, Florida 33710  
1-800-537-2790 (North America) • +32.11.27.04.17 (Europe) • (001)727 384-2323 (International)  
www.boviemed.com • info@boviemed.com

# MODULAR ERGONOMIC GRIP

**THE NEW STANDARD OF COMFORT,  
MODULARITY, SAFETY & ECONOMICS  
IN LAPAROSCOPIC SURGERY**

*Innovation  
At Its Best!*



**UNMATCHED WEALTH OF STYLES,  
DIAMETERS & USAGE MODES**



# A FRESH NEW LEVEL OF...

## ● ERGONOMICS

Ergonomically shaped, spring loaded Handle with ON-Demand Ratchet, reduces user fatigue.

## ● ECONOMICS

Single Handle serves 5mm working Diameters, Reusable or Disposable with different type jaws and scissors. Cost & Inventory are reduced.

## ● VERSATILITY

High Quality Reusable Long Lasting TIP-ON-TUBE™ Inserts with Brush-Thru Tube or Disposable Cartridges.

## ● SAFETY

Single Handle serves 5mm working Diameters, Metallic Handle is well insulated for safe, enduring, quality feel.

## ● ECO-FRIENDLY

Waste Reduction with Long Lasting Reusable Handle & Reusable Inserts.



## ADDITIONAL FEATURES THAT CAN SAVE YOU TIME & MONEY

- Accepts insert diameters from 2mm - 12mm (both "hot" and inert tips making it the ideal handle for laparoscopy, or any minimally invasive procedure)
- Tip inserts "click" securely into place
- Metallic handles provide long-lasting durability and solid feel while body provides electrical insulation
- The handle's ratchet feature optimize procedures requiring graspers or dissectors.
- One of a kind spring load automatically returns jaws to open position further reducing hand fatigue
- Easy to clean and re-sterilize
- 360 degree rotation for easy access to all anatomy
- Precisely balanced to accommodate any size hand with comfort, strength and agility

# 2 TIP OPTIONS PROVIDE CONVENIENCE, SAFETY & VALUE

MEG™ laparoscopy tips come as convenient single use cartridges or with an easy to clean, reusable tube and insert. Both styles are competitively priced and are made from the highest quality materials and workmanship. The MEG™ handle allows you to easily switch between the two options so that you can "mix and match" the options, even in the same procedure.

## ● DISPOSABLE CARTRIDGES

The totally disposable cartridges "click" directly into the MEG™ handle and fully rotate.



## ● TIP-ON-TUBE™

The reusable tube also "clicks" directly into the MEG™ handle and fully rotates. All tips simply "click" into the tube and are easily released. The outside is insulated and the inside is hollow for easy cleaning.