

ELECTROSURGERY FOR ANIMAL HEALTH



A1250S-V



A2350-V



A3350-V



BENEFITS	A1250S-V	A2350-V	A3350-V
Power (in watts)	120	200	300
Cut modes available	1	3	3
Blend modes (# of levels) available	1	4	4
Coagulation modes available	2	3	3
Bipolar modes available	1	3	4
Active cable connection	✓	✓	✓
Compatible with cost-saving reusable accessories	✓	✓	✓
Simplified transition from cut to coag on the autoclavable and reusable handpiece	✓	✓	✓
Fast Digital Feedback System (FDFS) - provides constant (non-fluxuating) energy through skin, fat or muscle	✓	✓	✓
Gentle Coag mode produces faster post-op wound healing		✓	✓
Time-saving presets for quick and easy setup		✓	✓
Laparoscopic mode reduces risk of alternate site burns with safety-enhanced lower voltage		✓	✓
Dual Pencil Port		✓	✓
Separate footcontrol for bipolar activation		✓	✓
Coag to bipolar transition in sterile field without ESU contact		✓	✓
Bovie Bipolar, a pulsating energy, can be used for fusing tissue together			✓
Dedicated Bipolar Presets			✓
Auto Bipolar - once resistance level is set, the energy level stops automatically. No bipolar footswitch is required with this mode.			✓
Warranty	4 years	4 years	4 years

Smoke Shark™ II

Animal health professionals can be exposed to surgical smoke proven to carry viruses, CO, and other chemical and biological hazards.









The Smoke Shark II minimizes exposure to these toxic smoke plumes during electrosurgery, laser, and harmonic scalpel procedures.



Bovie ANIMAL HEALTH

Ph +1 727 384 2323
 Sales@BovieMed.com
 www.BovieVet.com

ANIMAL HEALTH GENERATOR PACKAGES INCLUDE:

			
			
<p>A1250S-V</p> <ul style="list-style-type: none"> • A1250S-V Electrocauter Generator • One Reusable Pencil (ESPR2) • Five Non-Sterile Reusable Electrodes (Two ES01R Blades, Two A833 Needles, & One ES20R Ball) • One Metal Plate (A1204P) • Grounding Cord (A1254C) • Power Cord* • User's Guide* <p>SUGGESTED ACCESSORIES*</p> <ul style="list-style-type: none"> • Mobile Stand (ESMS-C) • Smoke Shark™ II Smoke Evacuator (SE02) • Monopolar Footswitch (BV-1253B) • Bipolar Footswitch (BV-1254B) • Footswitch Adapter (A1255A) • Bipolar Forceps Cable (A827V) • Bipolar Forceps • Suction Coagulators • Reusable SuperCut™ Needle Electrodes • Reusable Dermal Loops 	<p>A1250S-VS</p> <ul style="list-style-type: none"> • A1250S-V Electrocauter Generator • One Sterile Reusable Pencil (ESPR2) • Five Non-Sterile Reusable Electrodes (Two ES01R Blades, Two A833 Needles, & One ES20R Ball) • One Metal Plate (A1204P) • Grounding Cord (A1254C) • Mobile Stand (ESMS-C) • Smoke Shark II Smoke Evacuator (SE02) • 35-Hour Filter (SF35) • Pneumatic Footswitch (FSSEP)* • Two Power Cords* • User's Guides* 	<p>A2350-V</p> <ul style="list-style-type: none"> • A2350-V Electrocauter Generator • One Reusable Pencil (ESPR2) • One Metal Plate (A1204P) • Grounding Cord (A1254C) • Power Cord* • User's Guide* <p>SUGGESTED ACCESSORIES*</p> <ul style="list-style-type: none"> • Mobile Cart Stand (BV-IDS-CS) • Smoke Shark II Smoke Evacuator (SE02) • Monopolar Footswitch (BV-1253B) • Bipolar Footswitch (BV-1254B) • Footswitch Adapter (A1255A) • Bipolar Cord (A827V) • Bipolar Forceps • Suction Coagulators • Reusable SuperCut Needle Electrodes • Reusable Dermal Loops 	<p>A3350-V</p> <ul style="list-style-type: none"> • A3350-V Electrocauter Generator • One Reusable Pencil (ESPR2) • One Metal Plate (A1204P) • Grounding Cord (A1254C) • Power Cord* • User's Guide* <p>SUGGESTED ACCESSORIES*</p> <ul style="list-style-type: none"> • Mobile Cart Stand (BV-IDS-CS) • Smoke Shark II Smoke Evacuator (SE02) • Monopolar Footswitch (BV-1253B) • Bipolar Footswitch (BV-1254B) • Footswitch Adapter (A1255A) • Bipolar Cord (A827V) • Bipolar Forceps • Suction Coagulators • Laparoscopic Electrodes



Mobile Generator Stand
(BV-IDS-CS)



Monopolar
Footswitch
(BV-1253B)



Bipolar
Footswitch
(BV-1254B)

*Not pictured

Rev. 5