

## MI-300 MI-150

### Worth closer examination

#### LED Examination Lights



Presenting the best in high-value, space-saving, color-accurate LED exam lights. The MI-300 offers two brightness settings and its mobile base allows it to be moved anywhere easily. Both exam lights provide brilliant, color accurate illuminating and clarity.

#### SPECIFICATIONS

##### Light Intensity

MI-300.....	30,000 or 15,000 lux at 24" (60cm) 50,000 lux or 25,000 lux at 18" (45.7cm)
MI-150.....	15,000 lux at 24" (60cm) 25,000 lux at 18" (45.7cm)

Color Temperature.....4,300° K

CRI.....97

LED Type/Life.....CoB / 50,000 hours

Operating Voltage.....100VAC-240VAC/50-60Hz

Light Head Input Voltage.....24VDC


Unit Height.....MI-300: 66.0" (1,625mm)

MI-150: 64.0" (1,676mm)

Warranty.....3 years (limited)



Dual Luminance Settings



Mobile Base MI-300  
Stationary Base MI-150



Brilliant, Color-Accurate Illumination

Visit [BovieMedical.com](http://BovieMedical.com) for more information

#### LED Examination Lights

CODE	DESCRIPTION
XLDEM-150	MI-150 Exam Light Connectible Goose Neck w/Fixed Base (061223)
XLDEM-300	MI-300 Exam Light Flexible Goose Neck w/Rolling Base (061123)